

Transformation Protocol Using Heat Shock

Mary Teruel, 9/7/00

- 1) Take out competent E.coli cells from -80°C freezer.
- 2) Turn on water bath to 42°C .
- 3) Put 30-50 ul of competent cells in a 1.5 ml tube (Eppendorf or similar) – you may need more or less cells, depending how competent they are
- 4) Keep tubes on ice.
- 5) Add 50 ng of circular DNA into E.coli cells. Incubate on ice for 10 min to thaw competent cells.
- 6) Put tube(s) with DNA and E.coli into water bath at 42°C for 45 seconds.
- 7) Put tubes back on ice for 2 minutes to reduce damage to the E.coli cells.
- 8) Add 1 ml of LB (with no antibiotic added). Incubate tubes for 1 hour at 37°C .
(Can incubate tubes for 30 minutes, unless are trying to grow DNA for ligation which is more sensitive. For ligation DNA, leave tubes for 1 hour)
- 9) Spread about 100 ul of the resulting culture on LB plates (with the appropriate antibiotic added – usually Ampicillin or Kanamycin). Grow overnight.
- 10) Pick the colonies about 12-16 hours later.

-